

SUPINE POSITION - INSTRUCTIONS FLEXION ABDUCTION



- If the Duorole does not offer sufficient support for non-healing decubitus between the knees, near the tailbone, or seat bones in combination with leg spasms.

- CVA (stroke), joint disorders, MS (Multiple Sclerosis), and paraplegia.
- Back pain, muscular dystrophy, and pelvic instability.

Gently move the patient with a smooth hip movement on the Flexion Abduction or use the passive hoist. Make sure the top part (cover) is empty.



Slide the cover between the knees and insert the filling (Luna Pearl).



Use your hand to smooth the pearls between the knees.



The Flexion Abduction does not need to be removed for a hygienic cleaning of the groin or tailbone area. Simply pushing it back gives room for cleaning.

The patient is now in a comfortable lying position with the Flexion Abduction.

4



Make sure the pearls are divided evenly. Use your hand to smooth the pearls between the knees.



When the patient has long thighs you can place a Role under the Flexion Abduction.



When the patient has long lower legs, you can place a Luna or Crescent at the calves.