


SUPINE POSITION - COMPLETE SUPPORT (WITH LATERAL POSITION)

Distribute the pillow stuffing to ensure that it is divided evenly across the entire pillow. Use the **Comfort XL** for adults and the **Comfort** for children



Role and Minima can possibly be replaced by Duorole or Flexion Abduction

Leave one hand's breadth of space at end of Comfort XL after patient is positioned. This ensures that the tailbone is free



- High muscle tension
- Shoulder, spine, seat bone and/or heel decubitus
- Knee overstretching
 - Foot spasm
 - Comfort




Leave hand breadth's of space for tailbone between both ends of Comfort XL

Slide away one side of Comfort XL from underneath patient for 30-degree lateral position. Put Minima under knee or leave Role in place. If need be, place Luna Pearl between the knees



Fold line



1

No force (do not forcefully push or pull on the patient).

2

Always leave a hand's breadth of space.

3

Fill all of the spaces where the patient's body does not make contact with the pillow.

4

Smooth the area where the patient touches the pillow to relieve any air pressure.

5

Provide opposite support to the patient in a lateral position.