


SEATING SUPPORT - PREVENT SLUMPING




Distribute the pillow stuffing to ensure that it is divided evenly across the entire pillow. Use the **Comfort XL** for adults and the **Comfort** for children

- Neck problems (Crescent)
- Slumping in bed (Comfort XL)
- Pain in tailbone and seating bones when sitting
- Ability to remain seated longer and better
 - Assistance for mobilisation
- Support when falling over
 - Vertigo
 - Comfort

Leave three hand's breadths of space so that patient is lying half on Comfort XL and half in between



Fold line



- 1**
No force (do not forcefully push or pull on the patient).
- 2**
Always leave a hand's breadth of space.
- 3**
Fill all of the spaces where the patient's body does not make contact with the pillow.
- 4**
Smooth the area where the patient touches the pillow to relieve any air pressure.
- 5**
Provide opposite support to the patient in a lateral position.