

THE NUN



The Nun Is made of breathing material to make sure that the wound is not closed off.

The Nun cannot be physically attached to the shoe, but it does stay in place thanks to the perfect form.



THE NUN

The Nun is a useful tool that can be used to protect the bandage shoe by REBACARE® in the case of exuding wounds. The Nun can be placed deep into REBACARE®'s bandage shoes, and functions as an extra layer between bandages and the shoe.

The Nun does not touch the skin. If desirable, the Nun can be sterilised. Using a ballpoint pen, the name of the user can be written on the washing label. The Nun cannot be physically attached to the shoe, but it does stay in place thanks to the perfect form.

THE NUN

- Is made of breathing material to make sure that the wound is not closed off
- Can be cleaned using 90% alcohol
- Is resistant for 70 % ethanol
- Can be sterilised

- Can be used for both the right and the left foot
- Adapts to the temperature of its environment, which causes it to feel comfortable
- Is friendly to the skin
- Can be used hygienically. It is both anti house dust mite and anti-allergic, which makes it extra suitable for health care centres, nursing homes and hospitals
- Can be washed on 95 degrees Celcius
- Is provided with a CE-hallmark
- Can be use in the bandage shoes by REBACARE®

THE NUN IS SUITABLE FOR:

- Private use
- Nursing homes
- Hospitals

Description: The Nun

Item number: 7772

Size: One size