


THE STABILO AND BLANKET CUSHION

The size can be adjusted to the patient using Velcro.



The Stabilo is made of two blocks and can be placed around the ankle as well as the knee.

The Blanket cushion can be placed at the back or laterally on the mattress.



THE STABILO

REBACARE[®]'s Stabilo is meant for preventive and curative use:

- For patients with a high chance of a sustaining club feet (pes equinus) as a result of endo- and exorotation.
- For patients with a high chance of sustaining club feet as a result of a CVA, paraplegia, casts and geriatrics.

THE STABILO

- The Stabilo is made of two blocks and can be placed around the ankle as well as the knee
- The Stabilo makes the use of large leg rests and other mobility impairing measures unnecessary
- The size can be adjusted to the patient using Velcro
- The Stabilo can be used for the left and the right knee and ankle

Description: The Stabilo

Item number: 7773

Size: One size

THE BLANKET CUSHION

REBACARE[®]'s Blanket cushion is meant for preventive and curative use:

- For patients with a higher risk on decubitus of the toe and heel and for patients with a higher risk of club feet as a result of pressure, immobility and spasms.

The Blanket cushion is placed at the back or laterally on the mattress. It can be positioned correctly by placing pillows between the Blanket cushion and the bed rail.

THE STABILO AND THE BLANKET CUSHION:

- Can be cleaned hygienically
- Is alcohol resistant (90 %)
- Is easily cleaned and put away
- Cannot damage beds and mattresses
- Have easily washable covers (60 degrees Celcius)
- Have been provided with a CE-Hallmark

Description: The Blanket cushion

Item number: 7774

Size: One size