

DYNAMIC POSITIONING ORTHOSIS (DPO)

THE DPO OF REBACARE[®] OFFER COMFORT

The DPO are meant for preventive and curative use. The DPO offer comfort for patients with a high risk of decubitus and muscle strain requiring a changing lying position.

Hygiene and incontinence have been taken into account: The seams have been welded and are provided with a breathing filter.

- It contributes to the well-being of the patients
- By positioning the patient correctly into the DPO, pressure on the body can be reduced
- Satisfies the EPUAP-Norm
- Is provided with a CE-hallmark
- Is used hygienically
- Can be cleaned with 70% ethanol or 90% alcohol
- Fits within the WIP guidelines
- Comfortable, soft cover
- The soft stretch covers are washable at 60° C and the 100% cotton covers at 90° C.

THE DPO ARE SUITABLE FOR:

- Private use
- Health care centres
- Nursing homes
- Hospitals
- Intensive care medicine



THE MAXIMA (STANDARD SIZED PILLOW)

Experience support and ultimate comfort on the pillow produced by REBACARE[®].

Description: The Maxima
Size: 50 cm x 65 cm
Item number: 7777



THE MINIMA

The Minima is a small and multifunctional pillow that can be put between the knees, legs, hands and head. It can also be used as a headrest.

Description: The Minima
Size: 35 cm x 25 cm
Item number: 7771



THE ROLE

The Role can offer support for a lateral position, and can prevent a patient from rolling too far.

Description: The Role
Size: 95 cm x 35 cm
Item number: 7778



THE DUOROLE

The Duorole can be used to expose the heels, and can prevent the development of decubitus caused by spasms.

Description: The Duorole
Size: (2x) 95 cm x 35 cm
Item number: 7778



THE COMFORT AND COMFORT XL

The Comfort made by REBACARE® leaves an open space around the back of the head, the shoulders and the elbows. The Comfort XL creates additional open space around the spine and tail bone.

Description: The Comfort and Comfort XL
Size: L (100 cm) and XL (150 cm)
Item number: 77710 (L) and 77710 (XL)



THE FREE EAR® EARPILLOW

The innovative Free Ear® earpillow leave the space around the ear open. The Free Ear® earpillow is patented.

Description: Free Ear®
Size: 50 cm x 65 cm
Item number: 7776



THE SUPPORT

The Support simplifies repositioning and supports a position lying on the stomach.

Description: The Support
Size: 205 cm
Item number: 7779

The pillows have been welded completely, which causes the seams to be impenetrable for bacteria and fluids.

The pillows have been provided with a breathing filter.



Stretch covers for the DPO can be ordered separately by stating the item number + H. And covers 100% cotton (washable at 90 °C) under mention of + HK.