

# THE HOLWI ARMREST


*Reduces decubitus in the elbow by providing a large supporting surface of high quality foam.*



*Centre of rotation can be adjusted to joint of the elbow.*



*Can be used for both the left and the right arm.*



## THE HOLWI ARMREST MADE BY REBACARE®

The Holwi armrest is used for people suffering from hemiplegia, spasms and/or the loss of several functions of the arm. The Holwi armrest assures a proper pose. It prevents the development of pain in the shoulder, lower arm and/or the hand.

The Holwi armrest can be fixed onto every wheelchair. A great advantage of the armrest is that it does not put any pressure on the stomach. This is caused by the anatomic design as there are no rims at the side of the stomach.

### THE HOLWI ARMREST:

- Reduces of decubitus in the elbow by providing a large supporting surface of high quality foam
- Can be used for both the left and the right arm
- Has a centre of rotation that can be adjusted to the joint of the elbow
- Is steplessly adjustable for both the angle and hollowness

- Can adjust the rolling resistance to personal needs
- Is biomechanic
- Can be pushed away effortlessly by the user
- Makes independent transfers possible again
- Prevents a loose hanging arm next to the wheelchair by providing an open rim
- Does not put any pressure on the stomach because of its anatomic design
- Is skin friendly and can be used hygienically
- Is easily cleaned with a damp cloth

### THE HOLWI ARMREST IS SUITABLE FOR:

- Private use
- Health care centres
- Nursing homes
- Hospitals

**Description: Holwi armrest**

**Item number: 77713**

**Size: One size for all wheelchairs**