

## NUN



*Nun is made of breathing material to make sure that the wound is not closed off.*

*Nun cannot be physically attached to the shoe, but it does stay in place thanks to the perfect form.*



### NUN

Nun is a useful tool that can be used to protect the Orthopedic Bandage shoe / heel protector by REBACARE® in the case of exuding wounds. Nun can be placed deep into the Orthopedic Bandage shoe / heel protector, and functions as an extra layer between bandages and the shoe.

Nun does not touch the skin. If desirable, the Nun can be sterilised. Using a ballpoint pen, the name of the user can be written on the washing label. Nun cannot be physically attached to the shoe, but it does stay in place thanks to the perfect form.

### NUN

- Is made of breathing material to make sure that the wound is not closed off
- Disinfect with 70% alcohol and max. 1000 ppm Chlorine, Oxyfoam, detergent or OxyWipe
- Can be used for both the right and the left foot

- Can be sterilised
- Adapts to the temperature of its environment, which causes it to feel comfortable
- Is friendly to the skin
- Can be used hygienically. It is both anti house dust mite and anti-allergic, which makes it extra suitable for health care centres, nursing homes and hospitals
- Can be washed on 95 degrees Celcius
- Is provided with a CE-hallmark
- Can be use in the Orthopedic Bandage shoe / heel protector

### NUN IS SUITABLE FOR:

- Private use
- Nursing homes
- Hospitals

**Description:** Nun  
**Item number:** 7772  
**Size:** One size  
**CE:** NL-CA002-2011-21717

