


STABILO AND BLANKET CUSHION

The size can be adjusted to the patient using Velcro.



Stabilo is made of two blocks and can be placed around the ankle as well as the knee.

Blanket cushion can be placed at the back or laterally on the mattress.



STABILO

Stabilo is meant for *preventive* and *curative* use:

- For patients with a high chance of a sustaining club feet (pes equinus) as a result of endo- and exorotation.
- For patients with a high chance of sustaining club feet as a result of a CVA, paraplegia, casts and geriatrics.

STABILO

- Stabilo is made of two blocks and can be placed around the ankle as well as the knee
- Stabilo makes the use of large leg rests and other mobility impairing measures unnecessary
- The size can be adjusted to the patient using Velcro
- Stabilo can be used for the left and the right knee and ankle

Description: Stabilo

Item number: 7773

Size: One size

CE: NL-CA002-2013-29283



BLANKET CUSHION

Blanket cushion is meant for *preventive* and *curative* use:

- For patients with a higher risk on decubitus of the toe and heel and for patients with a higher risk of club feet as a result of pressure, immobility and spasms.

Blanket cushion is placed at the back or laterally on the mattress. It can be positioned correctly by placing pillows between the Blanket cushion and the bed rail.

STABILO AND BLANKET CUSHION:

- Can be cleaned hygienically
- Disinfect with 70% alcohol and max. 1000 ppm Chlorine, Oxyfoam, detergent or OxyWipe
- Cannot damage beds and mattresses
- Have been provided with a CE-Hallmark

Description: Blanket cushion

Item number: 7774

Size: One size

CE: NL-CA002-2013-29284

