

## DYNAMIC POSITIONING ORTHOSIS (DPO)

*The Dynamic Positioning Orthosis (DPO) from REBACARE® minimises pressure and reduces pain. The DPO is designed for preventive and curative use. It provides comfortable support for patients who have an increased risk of pressure injury or muscle strain and therefore require repositioning.*

### THE DPO FROM REBACARE®

- The seams are fully welded and impenetrable to moisture and bacteria.
- It helps prevent cross-infection.
- The permeable filter repels moisture while allowing air to pass.
- It helps correct positioning without separate supports.
- Has CE marking

- Correct positioning relieves pressure on the body.
- Complies with EPUAP standards
- Disinfect with 70% alcohol and max. 1000 ppm Chlorine, Oxyfoam, detergent or OxyWipe
- Meets Netherlands Working Group on Infection Prevention (WIP) guidelines
- Comfortable soft cover
- The soft stretch covers can be washed at 60 degrees Celsius, the 100% cotton covers at 95 degrees.

### THE DPO FROM REBACARE® CAN BE USED FOR:

Private use, Care homes, Nursing homes, Hospitals and Intensive Care Units.

Stretch covers (80% cotton and 20% polyamide) can be washed at 60 degrees Celsius. The stretch covers for the DPO can be ordered separately, stating the relevant article number + H, the 100% cotton covers (washable at 95 degrees Celsius) can be ordered stating the article number + HK.

*The permeable filter repels moisture while allowing air to pass.*

*The seams are fully welded and impenetrable to moisture and bacteria.*

## VARIOUS TYPES OF STUFFING

1

### SILICONISED HOLLOW FIBRE AND EPS (PEARLS)

*Very good quality, retains volume and shape. Is very springy. A combination of two types (2 and 3) of stuffing. Offers support and effectively reduces pressure.*



2

### SILICONISED HOLLOW FIBRE

*The fibre is springy, hollow and light. Feels comfortable and does not change shape. Because of its resilience, it is highly effective in absorbing pressure and offers good support.*



3

### EPS (PEARLS)

*These identical EPS pearls (1 mm in size) ensure optimal support of the patient. The combination of air and pearls serves to reduce pressure. A body print is made to create support that is perfectly moulded to the body.*





### MINIMA

A small and multi-purpose pillow. Can be used as a headrest, between the knees, or to provide opposite support in lateral position.

Description: Minima

Sizes: width 35 cm x height 25 cm

Article number: 77711

CE: NL-CA002-2014-33425



### MAXIMA

Experience support and supreme comfort with the Maxima pillow.

Description: Maxima

Sizes: width 70 cm x height 60 cm

Article number: 7777

CE: NL-CA002-2014-33425



### LUNA


Can be used as head pillow and is suitable for providing opposite support under the knee. Can be used in combination with the Role as abduction pillow between the legs. The Luna is a handy help for 15-degree prone position. It keeps the scrotum and chest off the mattress. Can also be used as ankle support in supine position.

Description: Luna

Sizes: width 70 cm x height 30 cm

Article number: 77717

CE: NL-CA002-2014-33425



### FREE EAR® EAR PILLOW

Leaves the ear free. By placing A on B (fix if necessary), the back of the head is kept off the mattress.

Description: Free Ear® Ear pillow

Sizes: width 65 cm x height 50 cm

Article number: 7776

CE: NL-CA002-2014-33425



### LUNA PEARL

(CVA) hemiplegia arm rest for in bed or chair. Also effective in preventing pes equinus in supine position. Convenient to use under the elbows when leaning heavily due to COPD, for instance. Suitable as lower arm support in case of shoulder fracture or cuff rupture.

Description: Luna Pearl

Sizes: width 70 cm x height 30 cm

Article number: 77721

CE: NL-CA002-2014-33425



### FREE EAR® TUBE PRONE PILLOW

Leaves the ear free and is suitable for prone position. Provides room for the Tube. Leaves the ear free and prevents ear pressure injury. Places the head perfectly flat for good posture. Free Ear® Tube Prone pillow has been patented.

Description: Free Ear® Tube Prone pillow

Sizes: width 65 cm x height 50 cm

Article number: 7776 T

CE: NL-CA002-2014-33425



### ROLE

The Role is suitable for under the knees (fowler) and can prevent a patient from rolling too far or back. In combination with the Luna, the Role can be used as a flexion-abduction pillow. The Role can be used under the elbows when leaning heavily due to COPD, for instance.

Description: Role

Sizes: width 95 cm x height 35 cm

Article number: 7778

CE: NL-CA002-2014-33425



### DUOROLE

The Duorole pillow can be used to keep the heels off the mattress. It prevents chafing pressure injury. Together with the fowler, it removes the pressure from the tail bone and seat bones. Improves muscle tone and helps break through spasms.

Description: Duorole

Sizes: (2x) width 95 cm x height 35 cm

Article number: 7778

CE: NL-CA002-2014-33425



### COMFORT AND COMFORT XL

The Comfort creates an open space around the back of the head, shoulders and elbows. The Comfort XL creates an additional open space around the spine and tail bone. Both can be used to provide opposite support to prevent slumping.

Description: Comfort and Comfort XL

Size L: width 100 cm x height 35 cm

Size XL: width 150 cm x height 35 cm

Article number: 77710(L) and 77710(XL)

CE: NL-CA002-2014-33425



### SUPPORT

Support provides a 30-degree lateral position and is used together with the Minima pillow. Also a good support for 15-degree prone position.

Description: Support

Sizes: width 205 cm x height 35 cm

Article number: 7779

CE: NL-CA002-2014-33425



### FORMA

For good support along the entire body. The body is not corrected, but accepted.

Description: Forma

Sizes: 200 x 35 cm

Article number: 77716

CE: NL-CA002-2014-33425



### TWISTER

Helps to keep the tail bone, spine and heels off the mattress. Convenient for use in Intensive Care units. The Twister is put under the patient with minimal twisting. The Twister is often used in combination with the Luna. In flat prone position, it keeps the scrotum off the mattress.

Description: Twister

Sizes: width 235 cm x height 16 cm

Article number: 77715

CE: NL-CA002-2014-33425



### CRESCENT

Can be used to support the neck and as forehead pillow in prone position. Can also be used as ankle support in supine position.

**Description:** Crescent

**Sizes:** width 50 cm x height 25 cm

**Article number:** 77718

**CE:** NL-CA002-2014-33425



### ORTHOPEDIC BANDAGE SHOE /HEEL PROTECTOR

Designed for preventive and curative use for patients suffering from or at risk of heel pressure injury.

- Unique stuffing composition that creates the ideal pressure distribution
- Keeps the heel, malleolus and the Achilles tendon off the mattress
- Prevents heel pressure injury
- Takes away the pain locally
- Safe anti-slip sole
- Stabilises the leg in combination with the Stabilo (prevents buckling)
- Complies with EPUAP standards
- Appropriate for the fall prevention policy
- Has CE marking
- Can be washed at 60° up to 95° Celsius

**Description:** Orthopedic Bandage shoe /heel protector

**Article number:** 7771

**Sizes:** width 30 cm x height 25 cm

**CE:** NL-CA002-2011-21717



### STABILO

Designed for preventive and curative use. Prevents buckling of patients with a heightened risk of endorotation and exorotation. Together with the Blanket Block, it helps prevent pes equinus. Suitable for patients with a heightened risk of pes equinus due to buckling (e.g. after CVA, spinal cord plaster and in geriatrics).

**Description:** Stabilo

**Article number:** 7773

**Sizes:** width 10 x height 12 x deep 23 cm per block

**CE:** NL-CA002-2013-29283



### BLANKET CUSHION

Designed for preventive and curative use. For patients with a heightened risk of toe and heel pressure injury. For patients with a heightened risk of pes equinus due to pressure, immobility and spasms.

- Disinfect with 70% alcohol and max. 1000 ppm Chlorine, Oxyfoam, detergent or OxyWipe
- Does not damage beds and mattresses
- Has CE marking

**Description:** Blanket Cushion

**Article number:** 7774

**Sizes:** width 68 x height 35 x deep 12 cm

**CE:** NL-CA002-2013-29284



**REBACARE GIVES A 2-YEAR PRODUCT WARRANTY ON THE DPO**  
**REBACARE GIVES A 6-MONTH PRODUCT WARRANTY ON THE ORTHOPEDIC BANDAGE SHOE /HEEL PROTECTOR**