


FREE EAR[®] EAR PILLOW AND FREE EAR[®] TUBE PRONE PILLOW

The seams are fully welded and impenetrable to moisture and bacteria. The permeable filter repels moisture while allowing air to pass.



The innovative Free Ear[®] Ear pillow leaves the ear free.

Free Ear[®] Ear pillow

Filled with siliconised hollow fibre. This is resilient and highly effective in absorbing pressure.

Free Ear[®] Tube prone pillow

Filled with EPS pearls. The combination of air and pearls serves to reduce pressure.

FREE EAR[®] EAR PILLOW

The innovative Free Ear[®] Ear pillow from REBACARE[®] is patented and leaves the ear free. Its use will help your organisation improve prevalence rates.

- Complies with EPUAP standards
- Has CE marking
- Hygienic in use
- Seams are fully welded and impenetrable to moisture and bacteria
- The permeable filter repels moisture while allowing air to pass
- Prevents moisture leaking into the product
- Disinfect with 70% alcohol and max. 1000 ppm Chlorine, Oxyfoam, detergent or OxyWipe
- Meets Netherlands Working Group on Infection Prevention (WIP) guidelines
- Comfortable soft cover
- The soft stretch covers can be washed at 60 degrees Celsius, the 100% cotton covers at 95 degrees

FREE EAR[®] EAR PILLOW IS SUITABLE FOR:

- Private use
- Care homes
- Nursing homes
- Hospitals
- Intensive Care units

INDICATIONS:

- Ear decubitus
- Chondrodermatitis nodularis helioides
- Oxygen tube over the ear
- Tube-feeding or mask over the face
- Hearing aid
- Tube prone pillow during ventilation

Description: **Free Ear[®] Ear pillow**

Sizes: 50 x 65 cm

Article number: 7776

Cover article number: 7776+H

CE: NL-CA002-2014-33425

Description: **Free Ear[®] Tube prone pillow**

Sizes: 50 x 65 cm

Article number: 7776 T

Cover article number: 7776+H

CE: NL-CA002-2014-33425

