

ORTHOPEDIC BANDAGE SHOE / HEEL PROTECTOR



Uniquely compiled stuffing which creates an ideal distribution of pressure.



Leaves open space around the heel, malleolus and the achilles tendon.

REBACARE®'S BANDAGE SHOE / HEEL PROTECTOR IS MEANT:

- For *preventive* and *curative* use for patients with a risk of getting heel pressure ulcers (decubitus).
- To *bridge the time* during which patients cannot wear their own shoes because of wounds, phlebology or lymphatic disorders.

REBACARE®'S BANDAGE SHOE / HEEL PROTECTOR CAN SAVE EXTRA COSTS BECAUSE:

- It improves the patients mobility, which lessens the duration of time the patient needs care.
- It enables fast healing of a wound, which lessens the duration of time the patient needs care.

REBACARE®'s bandage shoe / heel protector can gain your institute a good prevalence score.

REBACARE®'S BANDAGE SHOE / HEEL PROTECTOR

- Prevents decubitus on the heel

- Has a uniquely compiled stuffing which creates an ideal distribution of pressure
- Leaves open space around the heel, malleolus and the achilles tendon
- Takes away local pain
- Has a safe non-skid sole
- Stabilises the leg when used with the Stabilo
- Satisfies the EPUAP-norm
- Contributes to fall prevention measures
- Is provided with a CE hallmark
- Washable from 60 to a maximum of 95 degrees
- Can be used for the right and left foot

THE BANDAGE SHOE / HEEL PROTECTOR IS SUITABLE FOR

- Private use
- Intramural care

Description: bandage shoe / heel protector

Item number: 7771

Size: One size

CE: NL-CA002-2011-21717

